

# Veteran's Administration Catalog



[evologicsamerica.com](http://evologicsamerica.com)

4766 Research Drive

San Antonio, TX 78240

PHONE 866-982-1511

FAX 866-982-1533

FSS #36F79721D0009



**Because we care more.**



# PRODUCT NUMBERS

PRODUCT CODE	DESCRIPTION
<b>DBM PUTTY &amp; CRUSH MIX</b>	
<b>EHG-02C</b>	DBM CRUSH-MIX 2.5 cc
<b>EHG-05C</b>	DBM CRUSH-MIX 5 cc
<b>EHG-10C</b>	DBMCRUSH-MIX 10cc
<b>EHG-01P</b>	DBM PUTTY 1cc
<b>EHG-02P</b>	DBM PUTTY 2.5cc
<b>EHG-05P</b>	DBM PUTTY 5cc
<b>EHG-10P</b>	DBM PUTTY 10cc
<b>DBM CUBES &amp; STRIPS</b>	
<b>EV73010</b>	DBM CUBE 10x10x10mm
<b>EV73012</b>	DBM CUBE 12x12x12mm
<b>EV73014</b>	DBM CUBE 14x14x14mm
<b>EV50205</b>	DBM BONE SPONGE 50x20x5mm
<b>EVS 26197</b>	DBM Sponge, 26 x 19 x 7mm
<b>DBM FIBER BOATS</b>	
<b>EVB100205</b>	DBM FIBER BOAT 100x20x5 mm
<b>EVB502005</b>	DBM FIBER BOAT 50x20x5 mm
<b>CANCELLOUS CHIPS &amp; CUBES</b>	
<b>EV10015</b>	CANCELLOUS CRUSHED 1-4mm, 15cc
<b>EV10030</b>	CANCELLOUS CRUSHED 1-4mm, 30cc
<b>EV11015</b>	CANCELLOUS CRUSHED 4-10mm, 15cc
<b>EV11030</b>	CANCELLOUS CRUSHED 4-10mm, 30cc
<b>TENDONS</b>	
<b>AT652FF</b>	ANTERIOR OR POSTERIOR TIBIALIS TENDON
<b>AT290FF</b>	ACHILLES TENDON
<b>AT292FF</b>	BTB TENDON, HEMI
<b>AT393FF</b>	GRACILIS
<b>AT392FF</b>	SEMITENDINOSIS FROZEN GRAFT/ TENDON
<b>AT293FF</b>	PRE-SHAPED BTB TENDON
<b>AT294FF</b>	PERONEUS LONGUS
<b>AMNION PATCH</b>	
<b>EAP11</b>	EVOLOGICS AMNION PATCH 1 x 1cm
<b>EAP22</b>	EVOLOGICS AMNION PATCH 2x2 cm
<b>EAP23</b>	EVOLOGICS AMNION PATCH 2 x 3 cm
<b>EAP44</b>	EVOLOGICS AMNION PATCH 4 x 4 cm
<b>EAP46</b>	EVOLOGICS AMNION PATCH 4 x 6 cm
<b>EAP48</b>	EVOLOGICS AMNION PATCH 4 x 8 cm
<b>AMNION FLUID</b>	
<b>EAF 05</b>	.5ml AMNION FLUID
<b>EAF 10</b>	1.0ml AMNION FLUID
<b>EAF 20</b>	2.0ml AMNION FLUID
<b>EAF 05C</b>	.5ml AMNION FLUID, CRYOPRESERVED
<b>EAF 10C</b>	1.0ml AMNION FLUID, CRYOPRESERVED
<b>EAF 20C</b>	2.0ml AMNION FLUID, CRYOPRESERVED




# PRODUCTS

## EVOGRAFT™ DBM PUTTY




EVOGRAFT™ DBM PUTTY is formulated to resist irrigation with each lot tested to confirm osteoinductive potential. Evologics offers multiple sizes to address a wide array of procedures, including foot and ankle, trauma and extremities, knees, hips, craniomaxillofacial, spine and more. DBM can be used to augment autologous cancellous or corticocancellous grafts, and used in craniofacial surgery, spinal fusion, hand and foot surgery, joint reconstruction and fracture repair.

## EVOGRAFT™ DBM CRUNCH



EVOGRAFT™ DBM Crush is 100% demineralized bone matrix cortical ground mixed with mineralized cancellous chips. Evograft DBM Crush is formulated to resist irrigation with each lot tested to confirm osteoinductive potential. DBM can be used to augment autologous cancellous or corticocancellous grafts, and used in craniofacial surgery, spinal fusion, hand and foot surgery, joint reconstruction and fracture repair.

## EVOGRAFT™ DBM STRIP



EVOGRAFT™ DBM Sponge Strip consists of 100% demineralized bone matrix made from a single piece of cancellous bone, exposing essential growth factors to promote new bone formation. Upon compression, the strips will expand to its original state, filling the void and minimizing the space between the graft and host bone. DBM can be used to augment autologous cancellous or corticocancellous grafts, and used in craniofacial surgery, spinal fusion, hand and foot surgery, joint reconstruction and fracture repair.



# PRODUCTS

## EVOGRAFT™ DBM FIBER BOATS



EVOGRAFT™ demineralized cortical bone fibers that are entangled and shaped into sizes engineered to complement specific surgical applications. This unique process creates an interconnected graft material that contains BMPs and other growth factors necessary for the promotion of new bone formation. Evologics allografts are flexible upon hydration and each allograft has a sterility assurance level of  $10^{-6}$  via e-beam sterilization.

## EVOGRAFT™ CANCELLOUS CRUSHED CHIPS & CUBES



EVOGRAFT™ DBM Cancellous Crushed Chips are designed to fill bony voids or gaps, providing a scaffold for bone ingrowth with the patient's own bone. A ready-to-use alternative to autograft for bone regeneration, Evograft Cancellous Crushed Chips can be used to augment autologous cancellous or corticocancellous grafts, and used in craniofacial surgery, spinal fusion, hand and foot surgery, joint reconstruction and fracture repair.

## EVOGRAFT™ TENDONS



Evologics EvoGraft™ allograft tendons are sterilized using the exclusive PASCO<sub>2</sub>® process that leaves allograft tissue with more tensile strength and elastic properties than traditional methods like gamma irradiation. Using supercritical carbon dioxide, unwanted compounds are extracted while preserving the tensile and elastic strength of the tissue.




# PRODUCTS

## AMNION PATCHES



Evologics Amnion EvoPatch products are single- and dual-layer amniotic tissue membrane minimally manipulated to protect the collagen matrix and its natural properties. The Evologics Amnion Patch is commonly used to reduce scar tissue formation, reduce inflammation at the surgical site, healing enhancements and barrier protection. Evologics offers a single-layer and dual-layer patch. Both contain the same growth factors, but the single-layer graft can provide more comfort due to its soft and pliable design.

## LIQUID AMNION



Evologics Fluid Amnion and EvoFlow products are versatile, powerful tissue that supports the development of the fetus in utero. It is used for various types of surgical procedures, and can also be used for wound care and soft tissue injuries. Evologics is proud to partner with industry leaders that use stringent sterilization processes designed to maintain biocompatibility and integrity. Liquid Amnion can be used in wound management, chronic and non-healing dermal wounds, cutaneous wound care, reconstructive surgery, burn care, trauma and more.